





ENDODONTICS

NSK offers endodontic file and tips to achieve quick and effective treatments for root canal enlarging, irrigation, and filling.

ENDODONTICS




E4 • Application for root canal cleaning and removal of foreign particles.




NSK Satelec MODEL E4 ORDER CODE Z217428

E4D • Application for root canal cleaning and root canal enlargement.



NSK Satelec MODEL E4D ORDER CODE Z217528 • Diamond coated

E7 • Removal of filling and foreign material.



NSK Satelec MODEL E7 ORDER CODE Z217307 • Water not required

E7D • Enlargement of root canal wall.



NSK Satelec MODEL E7D ORDER CODE Z217317 • Diamond coated

E8 • Removal of filling and foreign material.



NSK Satelec MODEL E8 ORDER CODE Z217308 • Water not required

E8D • Enlargement of root canal wall.



NSK Satelec MODEL E8D ORDER CODE Z217318 • Diamond coated


*For tips E5, E6, E7, E7D, E8, E8D, E11 and E12, please use the attached E Tip Replacement Wrench.

*Cross section diagram is the image of a point 1 mm from the end. The illustrated tip diagram is four times bigger than the actual size.

ENDODONTICS




E5 • Lateral condensation.




NSK Satelec MODEL E5 ORDER CODE Z217305 • Water not required

E15D • For locating root canals and removal of calcifications.



NSK Satelec MODEL E15D ORDER CODE Z217319 • Diamond coated

E6 • Lateral condensation.



NSK Satelec MODEL E6 ORDER CODE Z217306 • Water not required

*For tips E5, E6, E7, E7D, E8, E8D, E11 and E12, please use the attached E Tip Replacement Wrench.

*Cross section diagram is the image of a point 1 mm from the end. The illustrated tip diagram is four times bigger than the actual size.

Satelec® is registered trademarks of Acteon SA, France. Refer to the power guide on pages 64-65.

ENDODONTICS



E11 • 120° Angle Holder for U files



- For U files (Shank dia : 0.8 mm)
- For root canal cleaning
- For anterior teeth
- Includes E Tip replacement wrench
- *File is optional

NSK Satelec[®] MODEL E11 ORDER CODE Z217031

E12 • 95° Angle Holder for U files



- For U files (Shank dia : 0.8 mm)
- For root canal cleaning
- For posterior teeth
- Includes E Tip replacement wrench
- *File is optional

NSK Satelec[®] MODEL E12 ORDER CODE Z217032

U Files


Length : 33 mm
Shank dia. : 0.8 mm

- ISO15 • Pack of 6
- ISO20 • Pack of 6
- ISO25 • Pack of 6
- ISO30 • Pack of 6
- ISO35 • Pack of 6


- MODEL U files 33 mm / #15 ORDER CODE Y900062
- MODEL U files 33 mm / #20 ORDER CODE Y900063
- MODEL U files 33 mm / #25 ORDER CODE Y900064
- MODEL U files 33 mm / #30 ORDER CODE Y900065
- MODEL U files 33 mm / #35 ORDER CODE Y900066

Always refer to this reference chart to ensure correct and accurate proper file length to avoid file breakage.

	#15	#20	#25	#30	#35
E11	22 mm	24.5 mm	22 mm	22 mm	22 mm
E12	22 mm	24.5 mm	24.5 mm	24.5 mm	22 mm




Endo Kit E11 For anterior teeth



NSK Satelec[®] MODEL Endo Kit E11 ORDER CODE Y900146

- E11 : 1 pc. 120° angle holder for anterior teeth
- E Tip replacement wrench : 1 pc.
- U files 33 mm #15, #20, #25, #30, #35 : Pack of 6 each size

Endo Kit E12 For posterior teeth



NSK Satelec[®] MODEL Endo Kit E12 ORDER CODE Y900147


- E12 : 1 pc. 95° angle holder for posterior teeth
- E Tip replacement wrench : 1 pc.
- U files 33 mm #15, #20, #25, #30, #35 : Pack of 6 each size

RETROGRADE ENDO

The microscopic diamond coated tip allows effective removal of dentin at the apex. 2 types of angled tips (90°, 70°) are available for anterior teeth and another 2 types of curved tips (right and left) are available for molar teeth.




E30RD • Diamond coated tip allows efficient cavity preparation.




NSK Satelec MODEL E30RD ORDER CODE Z217212 • Diamond coated • For posterior teeth (right angled)

E31D • Diamond coated tip allows efficient cavity preparation.




NSK Satelec MODEL E31D ORDER CODE Z217210 • Diamond coated • For anterior and posterior teeth (70°)

E30LD • Diamond coated tip allows efficient cavity preparation.




NSK Satelec MODEL E30LD ORDER CODE Z217213 • Diamond coated • For posterior teeth (left angled)


E32D • Diamond coated tip allows efficient cavity preparation.



NSK Satelec MODEL E32D ORDER CODE Z217211 • Diamond coated • For anterior teeth (90°)



Retro Kit



NSK Satelec MODEL Retro Kit ORDER CODE Y900149
 • E30RD, E30LD, E31D, E32D : 1 pc. each
 • Tip holder : 1 pc.