





















**VariaSurg** Ultrasonic Surgery Tips (TiN coating) **Bone Surgery**

 <p>MODEL: <b>SG1</b> order code: <b>Z305101</b></p> <ul style="list-style-type: none"> <li>• Dots mark 3, 6 and 9mm from top of tip • Five teeth: 0.5mm thick</li> </ul>	 <p>MODEL: <b>H-SG1</b> order code: <b>Z305151</b></p> <ul style="list-style-type: none"> <li>• Dots mark 3, 6 and 9mm from top of tip • Five teeth: 0.5mm thick</li> </ul>
 <p>MODEL: <b>SG1A</b> order code: <b>Z305138</b></p> <ul style="list-style-type: none"> <li>• Dots mark 3, 6, 9, 12 and 15mm from top of tip • Five teeth: 0.5mm thick</li> </ul>	 <p>MODEL: <b>SG2R</b> order code: <b>Z305102</b></p> <ul style="list-style-type: none"> <li>• Right curved tip* • Three teeth: 0.5mm thick</li> </ul>
 <p>MODEL: <b>SG2L</b> order code: <b>Z305103</b></p> <ul style="list-style-type: none"> <li>• Left curved tip* • Three teeth: 0.5mm thick</li> </ul>	 <p>MODEL: <b>SGB</b> order code: <b>Z305109</b></p> <ul style="list-style-type: none"> <li>• Dots mark 3, 6 and 9mm from top of tip • Three teeth: 0.5mm thick</li> </ul>
 <p>MODEL: <b>H-SGB</b> order code: <b>Z305155</b></p> <ul style="list-style-type: none"> <li>• Dots mark 3, 6 and 9mm from top of tip • Three teeth: 0.5mm thick</li> </ul>	 <p>MODEL: <b>SG8A</b> order code: <b>Z305139</b></p> <ul style="list-style-type: none"> <li>• Dots mark 3, 6, 9, 12 and 15mm from top of tip • Three teeth: 0.5mm thick</li> </ul>
 <p>MODEL: <b>H-SG8R</b> order code: <b>Z305156</b></p> <ul style="list-style-type: none"> <li>• Right curved tip* • Three teeth: 0.5mm thick</li> </ul>	 <p>MODEL: <b>H-SG8L</b> order code: <b>Z305157</b></p> <ul style="list-style-type: none"> <li>• Left curved tip* • Three teeth: 0.5mm thick</li> </ul>
 <p>MODEL: <b>SG14R</b> order code: <b>Z305122</b></p> <ul style="list-style-type: none"> <li>• Right curved tip* • Five teeth: 0.5mm thick</li> </ul>	 <p>MODEL: <b>SG14L</b> order code: <b>Z305123</b></p> <ul style="list-style-type: none"> <li>• Left curved tip* • Five teeth: 0.5mm thick</li> </ul>
 <p>MODEL: <b>SG19</b> order code: <b>Z305135</b></p> <ul style="list-style-type: none"> <li>• Dot mark 3mm from top of tip • Five teeth: 0.5mm thick</li> </ul>	 <p>MODEL: <b>SG30</b> order code: <b>Z305137</b></p> <ul style="list-style-type: none"> <li>• Sharp edge • 0.5mm thick</li> </ul>
 <p>MODEL: <b>SG59</b> order code: <b>Z305141</b></p> <ul style="list-style-type: none"> <li>• Dots mark 3, 6 and 9mm from top of tip • Three teeth: 0.5mm thick</li> </ul>	 <p>MODEL: <b>SG68</b> order code: <b>Z305143</b></p> <ul style="list-style-type: none"> <li>• Dots mark 2 and 4mm from top of tip • Three teeth: 0.5mm thick</li> </ul>

**VariaSurg** Ultrasonic Surgery Tips (TiN coating) **Scraper**

 <p>MODEL: <b>SG3</b> order code: <b>Z305104</b></p> <ul style="list-style-type: none"> <li>• Triple sided edge spatula tip</li> </ul>	 <p>MODEL: <b>SG4</b> order code: <b>Z305105</b></p> <ul style="list-style-type: none"> <li>• Edge spatula tip</li> </ul>
 <p>MODEL: <b>SG5</b> order code: <b>Z305106</b></p> <ul style="list-style-type: none"> <li>• Rounded edge spatula tip</li> </ul>	

**VariaSurg** Ultrasonic Surgery Tips (TiN coating) **Extraction**

 <p>MODEL: <b>SG17</b> order code: <b>Z305132</b></p> <ul style="list-style-type: none"> <li>• 0.7mm thick</li> </ul>	 <p>MODEL: <b>SG16R</b> order code: <b>Z305133</b></p> <ul style="list-style-type: none"> <li>• Right angled tip* • 0.7mm thick</li> </ul>
 <p>MODEL: <b>SG18L</b> order code: <b>Z305134</b></p> <ul style="list-style-type: none"> <li>• Left angled tip* • 0.7mm thick</li> </ul>	




**VariaSurg** Ultrasonic Surgery Tips **Socket Lift (Crestal Approach)**

 <p>MODEL: <b>SCL1</b> order code: <b>Z305101</b></p> <ul style="list-style-type: none"> <li>• Internal irrigation</li> </ul>	 <p>MODEL: <b>SCL1D</b> order code: <b>Z305151</b></p> <ul style="list-style-type: none"> <li>• Internal irrigation • Diamond coated</li> </ul>
 <p>MODEL: <b>SCL2</b> order code: <b>Z305138</b></p> <ul style="list-style-type: none"> <li>• Internal irrigation</li> </ul>	 <p>MODEL: <b>SCL2D</b> order code: <b>Z305138</b></p> <ul style="list-style-type: none"> <li>• Internal irrigation • Diamond coated</li> </ul>
 <p>MODEL: <b>SCL3</b> order code: <b>Z305102</b></p> <ul style="list-style-type: none"> <li>• Internal irrigation</li> </ul>	 <p>MODEL: <b>SCL3D</b> order code: <b>Z305103</b></p> <ul style="list-style-type: none"> <li>• Internal irrigation • Diamond coated</li> </ul>
 <p>MODEL: <b>SCL4</b> order code: <b>Z305109</b></p> <ul style="list-style-type: none"> <li>• Internal irrigation • Diamond coated</li> </ul>	 <p>MODEL: <b>SCL5</b> order code: <b>Z305155</b></p> <ul style="list-style-type: none"> <li>• Internal irrigation</li> </ul>
 <p>MODEL: <b>SCL5D</b> order code: <b>Z305139</b></p> <ul style="list-style-type: none"> <li>• Internal irrigation • Diamond coated</li> </ul>	

**VariaSurg** Ultrasonic Surgery Tips (Diamond coating) **Sinus Lift**

 <p>MODEL: <b>SG6D</b> order code: <b>Z305107</b></p>	 <p>MODEL: <b>SG7D</b> order code: <b>Z305108</b></p>
--	--

**VariaSurg** Ultrasonic Surgery Tips **Sinus Membrane Detachment**

 <p>MODEL: <b>SG9</b> order code: <b>Z305110</b></p> <ul style="list-style-type: none"> <li>• Flat circular convex elevator • Angled at 90°</li> </ul>	 <p>MODEL: <b>SG10</b> order code: <b>Z305111</b></p> <ul style="list-style-type: none"> <li>• Flat circular convex elevator • Angled at 135°</li> </ul>
 <p>MODEL: <b>SG11</b> order code: <b>Z305112</b></p> <ul style="list-style-type: none"> <li>• Cone compressor</li> </ul>	




**VariaSurg** Ultrasonic Surgery Tips (Diamond coating) **Implant Preparation**

 <p>MODEL: <b>SG15A</b> order code: <b>Z305124</b></p> <ul style="list-style-type: none"> <li>• Diameter of the tip end: 0.7mm</li> </ul>	 <p>MODEL: <b>SG15B</b> order code: <b>Z305125</b></p> <ul style="list-style-type: none"> <li>• Diameter of the tip end: 1.3mm</li> </ul>
 <p>MODEL: <b>SG15C</b> order code: <b>Z305126</b></p> <ul style="list-style-type: none"> <li>• Diameter of the tip end: 0.9mm</li> </ul>	 <p>MODEL: <b>SG15D</b> order code: <b>Z305127</b></p> <ul style="list-style-type: none"> <li>• Diameter of the tip end: 1.3mm</li> </ul>
 <p>MODEL: <b>SG16A</b> order code: <b>Z305128</b></p> <ul style="list-style-type: none"> <li>• Diameter: 2.0mm</li> </ul>	 <p>MODEL: <b>SG16B</b> order code: <b>Z305129</b></p> <ul style="list-style-type: none"> <li>• Diameter: 2.6mm</li> </ul>






**VariaSurg** Ultrasonic Tips **Scaling**

 <p>MODEL: <b>G-1S</b> order code: <b>Z305113</b></p>
--





**VariaSurg** Ultrasonic Tips **Perio (Root Planing)**

 <p>MODEL: <b>P20-S</b> order code: <b>Z305114</b></p> <ul style="list-style-type: none"> <li>• Straight tip</li> </ul>	 <p>MODEL: <b>P25R-S</b> order code: <b>Z305115</b></p> <ul style="list-style-type: none"> <li>• Right curved tip*</li> </ul>
 <p>MODEL: <b>P25L-S</b> order code: <b>Z305116</b></p> <ul style="list-style-type: none"> <li>• Left curved tip*</li> </ul>	

**VariaSurg** Ultrasonic Tips **Perio-Control** **Maintenance (V-Tip)**

 <p>MODEL: <b>V10-S</b> order code: <b>Z305117</b></p> <ul style="list-style-type: none"> <li>• Plastic Tip is not included • Includes E-Tip replacement wrench</li> </ul>	 <p>MODEL: <b>V-P10</b> order code: <b>Y900184</b></p> <ul style="list-style-type: none"> <li>• Pack of 3 • V10-S holder is not included</li> </ul>
 <p>MODEL: <b>V-P11R</b> order code: <b>Y1002165</b></p> <ul style="list-style-type: none"> <li>• Right curved tip*</li> </ul>	 <p>MODEL: <b>V-P11L</b> order code: <b>Y1002166</b></p> <ul style="list-style-type: none"> <li>• Left curved tip*</li> </ul>
 <p>MODEL: <b>V-P12</b> order code: <b>Y1002167</b></p> <ul style="list-style-type: none"> <li>• Pack of 3 • V10-S holder is not included</li> </ul>	

**VariaSurg** Ultrasonic Surgery Tips (Diamond coating) **Endodontics**

 <p>MODEL: <b>E30RD-S</b> order code: <b>Z305118</b></p> <ul style="list-style-type: none"> <li>• For posterior teeth (right angled)</li> </ul>	 <p>MODEL: <b>E30LD-S</b> order code: <b>Z305119</b></p> <ul style="list-style-type: none"> <li>• For posterior teeth (left angled)</li> </ul>
 <p>MODEL: <b>E31D-S</b> order code: <b>Z305120</b></p> <ul style="list-style-type: none"> <li>• For anterior and posterior teeth (70°)</li> </ul>	 <p>MODEL: <b>E32D-S</b> order code: <b>Z305121</b></p> <ul style="list-style-type: none"> <li>• For anterior teeth (90°)</li> </ul>

**VariaSurg** Ultrasonic Surgery Tips **Tip Kits**

 <p>MODEL: <b>Basic H-S Kit</b> order code: <b>Y1002775</b></p> <ul style="list-style-type: none"> <li>• Contains: • H-SG1, SG4, SG6D, SG9, SG10, SG11 • Tip holder</li> </ul>
 <p>MODEL: <b>Bone Cut Kit</b> order code: <b>Y900688</b></p> <ul style="list-style-type: none"> <li>• Contains: • SG1, SG3R, SG4, SG2L, SG6D • Tip holder</li> </ul>
 <p>MODEL: <b>Sinus Lift Kit</b> order code: <b>Y900689</b></p> <ul style="list-style-type: none"> <li>• Contains: • SG1, SG3, SG6D, SG9, SG10, SG11 • Tip holder</li> </ul>
 <p>MODEL: <b>Endo-S Kit</b> order code: <b>Y900691</b></p> <ul style="list-style-type: none"> <li>• Contains: • G1-S, E30RD-S, E30LD-S, E31D-S, E32D-S • Tip holder</li> </ul>
 <p>MODEL: <b>Implant Preparation Kit</b> order code: <b>Y900774</b></p> <ul style="list-style-type: none"> <li>• Contains: • SG15A, SG15B, SG15C, SG15D • Tip holder</li> </ul>
 <p>MODEL: <b>Socket LRF Kit</b> order code: <b>Y1002841</b></p> <ul style="list-style-type: none"> <li>• For Regular Size Implant</li> </ul>
 <p>MODEL: <b>Socket LRF Kit</b> order code: <b>Y1002842</b></p> <ul style="list-style-type: none"> <li>• For Wide Size Implant</li> </ul>

MODEL	PRODUCTS	ORDER CODE	MODEL	PRODUCTS	ORDER CODE	MODEL	PRODUCTS	ORDER CODE
E-Tip Replacement Wrench For V10-S		Z17139	Tip Wrench (E-Tip)		1000397	VS Tip Wrench (VS-E30RD-S/SG15)		Z305350
			Tip Holder		2001127			