




## Facebows

### Springbow Earpiece and Facia Facebows

The laboratory needs data to replicate a patient's functional patterns.

The facebow record provides the proper orientation of the maxillary model to the articulator. Without the facebow record, an articulator is ineffective in achieving proper occlusion.



### Springbow Earpiece Facebow (complete)

#### Features:

- Used to capture and record a patient's maxillary arch and its relationship to the condyles
- Quick and easy to use
- Allows indirect or direct mounting to Hanau™ and Denar® articulators (except Hanau™ Model 145 or Denar® Mark 300 Series)
- Automatically centers for outstanding accuracy and patient comfort, may be sterilized by steam autoclave

#### Comes complete with:

Springbow, transfer clamp assembly, anterior elevating assembly, mounting platform with built-in cast support, bitefork and biteplane

**Item: 010328-000**

### Springbow Direct Mounting Only (not shown)

#### Comes complete with:

Springbow, transfer clamp assembly, anterior elevating assembly and bitefork

**Item: 010351-000**

View utilization video  
online at:  
[www.whipmix.com](http://www.whipmix.com)