



Articulators: Wide-Vue

Wide-Vue Semi-Adjustable Articulator

- Lingual access is unobstructed. You'll appreciate its enhanced view when positioning tooth-to-tooth contacts on models and final prostheses.
- Closed tracking fossa and fixed upper member are ideal for fabricating prosthetic appliances.



**Customer Service:
1.800.626.5651**



Wide-Vue Articulator

Adjustments:

Progressive Side Shift: 0°-30°
Protrusive angle: -20°-+60°

Comes complete with:

Orbitale indicator, dual-end incisal pin, 10 disposable mounting plates, instruction manual, carrying case

Item: 010885-000



**Maxillo Facial Articulator
with additional 25mm (1") vertical height**

Adjustments:

Progressive Side Shift: 0°-30°
Protrusive angle: -20°-+60°

Comes complete with:

Orbitale indicator, dual-end incisal pin, 2 pairs of metal mounting plates, instruction manual

Item: 010882-000

Magnetic mounting plates optional
See page 29.

View utilization video online at:
www.whipmix.com

Wide-Vue Accessories



Adjustable Incisal Guide

May be moved anteriorly or posteriorly which minimizes having to adjust for the "x" dimension. Incline may be adjusted and locked from 0° to 40° protrusive. Adjustable lateral wings establish the degree of incisal freedom from 0° to 45°.

Item: 010731-000



Dual-End Incisal Pin

Has a chisel end for use with an adjustable incisal guide and a spherical end for use with a flat anterior guide. Both ends are calibrated for vertical adjustments.

Item: 010739-000



Adjustable Incisal Pin

Allows the vertical distance on the articulator to be altered within a range of 5mm. Remains on the exact rotational center of the incisal guide, permitting an increase or decrease of the guide inclination without altering the vertical dimension.

Item: 010670-000



Remounting Record Jig (Bonwill)

Used in preparation of a permanent plaster record of full or partial dentures on an articulator, maintaining original relationship to the condyles, and for "milling in" of finished dentures.

Item: 003102-000



Remounting Record Jig

Used for preparation of a permanent plaster record, accurately relocates the occlusal index for remounting the upper model after processing. This precludes taking a second facebow record and transfer to the articulator.

Item: 009544-000



Orbitale Indicator (Frankfort Horizontal Plane)

Used with procedures employing the infra-orbital notch as an anatomical third reference point during facebow application and transfer.

Item: 010046-000