



BEGO mould formers

- Eliminates fixing and grinding of the investment models when the BEGO combination duplicating flask is used
- Both mould formers can also be used with all other duplicating systems
- Bases for lifting technique, ideal for plotted CAD/CAM frames and partial denture frames made from light-curing wax

5

Product details

Availability	Contents	REF
Mould former, small, red	4 pieces	52390
Mould former, large, blue	4 pieces	52400
Base with funnel former set (2 × large / 2 × small)	4 pieces	52401
Silicon mould former incl. funnel former	1 piece	54877

