



HARVARD[®]

Harvard SuperSnap[®]

Vinylpolysiloxane impression material



Fast
as the
fastest!



EN

Impression material of a new Generation...

...that is Harvard SuperSnap®. This product combines all the advantages of our previous Premium-Sil products with new setting properties.

You are looking for a combination of a long working time and a short oral setting-time

- This allows a stress-free workflow, time savings and increases patient comfort.

You want to use two different types of setting times depending on the clinical situation.

- Available in 2 snap-set variants: **Regular Snap** and **Quick Snap**







You want an exact readability.

- Harvard SuperSnap® offers new high-contrast colors for exact readability.

Harvard SuperSnap® offers as usual:

- Excellent flow properties
- High elastic recovery and dimension stability
- Excellent legibility
- Super hydrophilicity for perfect results and excellent reproduction of details
- As usual: very berry!
- Suitable for the sandwich technique and the putty-wash technique
- Different viscosities: Putty, Putty Soft, Heavy, Mono, Medium, Light
- Optimal kneadable consistency for the putty products

The recommended combinations of Harvard SuperSnap®:

Technique	Phase 1	Phase 2
2 steps / 2 phases putty-wash technique (Material application successively)		
1 step / 2 phases sandwich technique (Material application simultaneously)	 <p>or</p> 	
1 step / 1 phase Monophasetechnique		

Numbers that speak for themselves!

Product	Working time 23°C / 73°F (min: sec.)		Oral setting time (min: sec.)	
	Quick Snap	Regular Snap	Quick Snap	Regular Snap
Harvard SuperSnap® Light	1:10	1:45	1:15	1:45
Harvard SuperSnap® Medium	1:10	1:45	1:15	1:45
Harvard SuperSnap® Mono	1:10	1:45	1:15	1:45
Harvard SuperSnap® Heavy	1:10	1:45	1:15	1:45
Harvard SuperSnap® Heavy MaxiMix	1:10	1:45	1:15	1:45
Harvard SuperSnap® Putty Soft MaxiMix	-	1:45	-	1:45
Harvard SuperSnap® Putty Handmix	1:10	1:45	1:15	1:45
Harvard SuperSnap® Putty Soft Handmix	1:10	1:45	1:15	1:45
	Fast Set	Normal Set	Fast Set	Normal Set
Harvard PremiumSil (Light, Medium, Mono, Heavy)	1:10	1:45	1:30	2:30
Harvard PremiumSil (Putty, Putty Soft)	1:00	1:20	1:45	2:30
Harvard PremiumSil (Putty Soft MaxiMix)	-	1:20	-	2:30
	Super Quick	Regular Set	Super Quick	Regular Set
3M Imprint 4 Ultra Light	1:15	-	1:15	-
3M Imprint 4 Light	1:15	2:00	1:15	2:00
3M Imprint 4 Regular	1:15	2:00	1:15	2:00
3M Imprint 4 Penta Heavy	1:15	2:00	1:15	2:00
3M Imprint 4 Penta Putty	-	1:30	-	2:30



Quick Snap

Order information

Harvard SuperSnap®

Regular Snap

2 automix cartridges ea. 50 ml, 12 mixing tips 1:1, 4 IntraTips yellow

Light orange 7083801
Medium green 7083811

2 automix cartridges ea. 50 ml, 12 mixing tips 1:1

Mono blue 7083821
Heavy black-grey 7083831

2 jars ea. 600 g, 2 spoons

Putty black 7083841
Putty Soft black 7083851

380 ml maximix cartridge, 10 dynamic mixers, 1 bayonet ring

Putty Soft black 7083843
Heavy black 7083833

Quick Snap

2 automix cartridges ea. 50 ml, 12 mixing tips 1:1, 4 IntraTips yellow

Light orange 7083802
Medium green 7083812

2 automix cartridges ea. 50 ml, 12 mixing tips 1:1

Mono blue 7083822
Heavy black-grey 7083832

2 jars ea. 600 g, 2 spoons

Putty black 7083842
Putty Soft black 7083852

380 ml maximix cartridge, 10 dynamic mixers, 1 bayonet ring

Heavy black-grey 7083834

Mixing Tips und Zubehör

Harvard Auto-T 1:1 O-Yellow 7098000
Refill bag with 50 Mischkanülen for
Light, Medium – Regular / Quick

Harvard Auto-T 1:1 O-Green 7098010
Refill bag with 50 Mischkanülen for
Heavy – Regular / Quick

Harvard Auto-T 1:1 O-Pink 7098020
Refill bag with 50 Mischkanülen for
Mono – Regular / Quick

Harvard IntraTips yellow 7083619
Refill bag with 50 intra tips for
Light, Medium – Regular / Quick

Harvard Maxi 5:1 Dynamic 7091400
Refill bag with 50 dynamic mixers and
2 bayonet rings for Putty Soft MaxiMix
Heavy – Regular / Quick MaxiMix

Harvard Dispenser Automix 1:1 7095100

Harvard TraySive 7083753
10 ml bottle



Liability is excluded for all printing errors and omissions. Before using our Harvard products, the respective directions for use should be noticed in every case. All measurements are internal measurements of Harvard Dental International.

Harvard Distribution Partner.



HARVARD®

Marke und Qualität seit 1892

Harvard Dental International GmbH
Margaretenstr. 2 - 4, 15366 Hoppegarten, Germany
Phone: + 49 (0) 30/99 28 978-0
Fax: + 49 (0) 30/99 28 978-19
info@harvard-dental-international.de