



BegoSol®

Mixing liquid for BEGO investment materials

- Depending on the alloy and the field of application, the required mixing ratio can be created for these liquids using distilled or demineralised water
- The higher the concentration of the mixing liquid, the greater the expansion of the investment material

5

Product details

| Availability | Contents | REF |
|---|--------------|-------|
| BegoSol®* Mixing liquid for Wirovest®plus, WiroFast, Wiroplus® S, Wirovest®, Bellavest® T and Bellasun | 1 l bottle | 51090 |
| BegoSol® | 5 l Kanister | 51091 |
| BegoSol® HE** Special-Mixing liquid for Bellavest® SH, Bellavest® DR, Bellavest® T, VarseoVest P | 1 l bottle | 51095 |
| BegoSol® HE | 5 l Kanister | 51096 |
| BegoSol® K** Special-Mixing liquid for WiroFine, BellaStar XL, VarseoVest P plus | 1 l bottle | 51120 |
| BegoSol® K | 5 l Kanister | 51121 |
| BegoSol® CC Special-Mixing liquid for VarseoVest C&B | 1 l bottle | 54907 |
| BegoSol® CC Special-Mixing liquid for VarseoVest C&B | 5 l Kanister | 54908 |
| Accessories | | |
| Universal measuring cup 100 ml | 1 piece | 14607 |

* Anti-freeze optimization up to -10 °C

** Is sensitive to freezing.

Images and illustrations are examples. Colors, symbols, designs, and information on the depicted labels and/or packaging may differ from reality.