



SGL70M / 81.9 mm, 280 g*

SGL50M / 98.1 mm, 322 g*

Well-balanced, with a Compact and Lightweight Micromotor

Streamlining efforts reduced the size by 16.2 mm and the weight by 42 g to improve balance during use the Surgic Pro motor, greatly reducing any strain on clinicians.

LED Optics for Safer, More Accurate Treatment



NSK LEDs generate natural daylight-quality light to illuminate the treatment area, enabling more precise surgery and shortened operation times. The lights increase safety because they do not overheat and are long-lasting.



Solid titanium body with Cord 2 m

