

**Customer Service:
1.800.626.5651**



Earpiece Facebow

Features:

- Suspension from the external auditory meatus makes it a one operator bow
- Meets the precise requirements for fixed restorative techniques
- Available in several styles to comply with various methods of facebow registration and transfer and is adaptable to most makes of articulators
- Scales are tilted forward for ease of reading and adjustment
- The Bitefork/Biteplane clamp has an overhead integral lock—no separate wrenches are required. Finger pressure achieves locking without introducing torque into the bow.

Comes complete with:

Orbitale pointer and bitefork

For all current Hanau™ articulators

Item: 008824-000



Facia Facebow

Features:

- Includes a pointer for alignment to the orbitale as a third point of reference
- Equipped with two clamps; one to retain the bitefork and one to retain the orbitale pointer

Comes complete with:

Orbitale pointer and bitefork

For all current Hanau™ articulators

Item: 005705-000



Facebows Accessories



Transfer Clamp Assembly

Item: 010356-000



Bitefork (edentulous use)

Item: 005726-000



Biteplane (dentulous use)

Item: 005727-000



Mounting Platform

Has a built-in cast support. Attaches to the articulator to transfer patient's maxillary record captured with bitefork assembly on the Springbow.

Item: 010342-000



Maxillary Cast Support

The adjustable "T" supports the bitefork and maxillary cast when mounting the facebow record on the articulator. Works with both magnetic and traditional plates.

Item: 003401-000

Snaps on magnetic system!