

Orthodontic Plaster

Super-White Dental Plaster



Super White

Developed specifically for orthodontic use, Orthodontic Plaster is harder and whiter than ordinary plaster. It has a long working time and trims easily. To provide an accurate and dense cast,

finest results are obtained when mixed under vacuum. To produce a glossy surface, finished models may be treated with Model Glow model soap.

ISO TYPE 2

Physical Properties

Water/Powder Ratio	37 mL/100 g
Working Time	7–9 minutes
Setting Time	14 minutes
Setting Expansion	0.20%
Compressive Strength, Wet (1 hour)	3,500 psi (25 MPa)
Compressive Strength, Dry (48 hours)	6,500 psi (45 MPa)

Item No.

Description

Super-White

03077	9 kg (20 lb.) Carton
31836	22 kg (50 lb.) Carton