



My vis-à-vis

«We and the same?! Even our flosses differ. One loves mint, the other doesn't. One uses the fluffy, the other the compact. We go shopping together anyway!»

Lauren & Charlotte,
(Identical Twins)



paro[®] tapes and flosses

paro[®] brush'n floss

The *paro*[®] brush'n floss is ideal for cleaning braces, bridges and wide gaps between teeth. Its unique components - a stiffened-end dental floss threader and spongy floss - all work together for maximum benefits.

paro[®] glide-tape

The *paro*[®] glide-tape made from PTFE is a proven oral hygiene aid. Used daily, it can remove plaque in tight interdental spaces which are unreachable by other means. And it is easy on the hands; it neither cuts into fingers nor cuts off circulation.

paro[®] riser-floss

paro[®] riser-floss is a medicinal dental floss. It is minted, waxed and fluoridated. It is thin before applying and gets fluffy while used.

paro[®] classic-floss

paro[®] classic-floss is a thin, waxed, minted and fluoridated dental floss. As a result, it slides easily between the teeth and therefore recommended for narrow interdental spaces.

paro[®] floss holder

paro[®] floss holder makes its application to a mere child's play for everyone. It is a uniquely vigorous tool for perfect cleaning effect. All interdental spaces can be reached without any difficulties. High hygiene standard by not touching the floss. Highest efficiency with small consumption. To use with *paro*[®] classic-floss or *paro*[®] riser-floss.

1760



1764



1762



1766



1768



1767



1765



1763



Art.-No. | paro® tapes and flosses

1760	paro® brush'n floss – 20x15 cm
1762	paro® glide-tape – PTFE-tape, 20 m
1768	paro® glide-tape – PTFE-tape, 200 m
1764	paro® riser-floss – waxed, mint, fluoridated, 50 m
1763	paro® travel-floss – waxed, mint, 5 m
1766	paro® classic-floss – waxed, mint, fluoridated, 50 m
1767	paro® classic-floss – waxed, 200 m
1765	paro® floss holder