



WiroGel® M

Environmentally friendly duplicating gel based on agar-agar hydrocolloid for producing models using investment material, plaster and the acrylic casting technique

- For universal use: For all phosphate-bonded investment materials and type 4 plasters as well as the acrylic casting technique
- Highly accurate impression-taking; smooth model surfaces guarantee reliability in use and enable work results which meet the highest demands in terms of precision
- 15 melting cycles mean a very good cost-benefit ratio – suitable for melting in a microwave without compromising quality or precision
- Duplicating with WiroGel® M is more than five times less expensive than with silicone, taking into account the costs for the duplicating unit
- Color geared to contrast optimisation, thus ensuring optimal process reliability

Product details

Availability	Contents	REF
WiroGel® M	6 kg tub	54351
WiroGel® M	10 kg tub	54354
Accessories		
Combi duplicating flask	1 set	52090